# Operation Manual Shodex KW402.5-4F

J-22

(Please read this manual carefully before using the column to keep its good performance and life.)

## 1. Introduction

Shodex KW402.5-4F is a size exclusion chromatography column developed for separation of hydrophilic substances with low molecular weight, such as small proteins and peptides.

## 2. Instructions in handling <Important>

### Caution!

\* Take notice of keeping instructions about the solvents and the reagents used with the column not to occur problems related to losing your health or leaking.

### Attention!

\* Use the column within the regular range of flow rate, pressure and temperature. There is a danger of deteriorating the performance when it is handled beyond the permissible range even for a short time. See the clause "Usable conditions" about the permissible range.

## 3. Specifications

Column size	: 4.6 mm I.D. × 300 mm L.
Column material (inside)	: SUS 316
Packing material	: Porous silica particle (3 μ m) bound with hydrophilic polymer
In-column solvent (initial)	: Water
Exclusion limit M.W.	: 60,000 (pullulan)
Number of theoretical plates	: ≥35,000 per column

## 4. Usable conditions

Flow rate	: ≤0.5 ml/min, normally used under 0.35ml/min. * It should be 0.15 ml/min or less for replacing the in-column solvent.
Pressure	: ≤10MPa per column
Temperature	: 5~45 °C
pH	: 3.0~7.5
Eluent	: An aqueous solution of salt or a buffer solution is usually used. <Representative salts> sodium chloride, sodium nitrate, sodium sulfate, potassium sulfate, ammonium sulfate. Acceptable concentration of total salts is 0.1M~0.3M. <Representative buffers> phosphate buffer, tris buffer, acetate buffer. * Methanol or acetonitrile can be added up to 100 vol%.

### Attention!

- 1) Do not remove the end fittings of the column under any circumstances.
- 2) Do not make a strong impact on the column: such as hitting or dropping on the floor.
- 3) Replace the solvent in the chromatograph with the eluent to be used before connecting the column.
- 4) Connect the column so that the flow direction corresponds to the arrow mark on the tag.
- 5) When the column is not used for a month or more, replace the in-column solvent with water, close each end with a stopper, and store it at room temperature.
- 6) Filtrate the sample with a disposable filter (0.45 μ m) to prevent deterioration by adsorbing insoluble matters.
- 7) When methanol or acetonitrile is added to an aqueous solution of salt or a buffer solution, pay attention not to separate any salt after mixing.

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